



NRV OrthoTech Pvt Ltd
Perfection At It's Finest...



Orthopedic Implants & Instruments Company

ACL PCL Reconstruction

Interlocking Nails

Nails / Wires / Pins

DHS / DCS Plates / Hip Prosthesis

Cannulated Screws

Spine Implants & Instruments

Maxillofacial Implant

Mini Plates - Standard

Small Plates - Standard

Large Plates - Standard

Mini Locking Plates

Small Locking Plates

Large Locking Plates

External Fixators / Tools

Innovative Reliable Quality
Excellent Performance



Distal Humeral Locking Plate, 2.7mm (L&R)

SSL14 / TT38

Type : Dorsolateral - with Lateral Support
Size : 3,5,7,9 & 14 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



Distal Humeral Locking Plate, 2.7mm (L&R)

SSL15 / TT39

Type : Dorsolateral - without Lateral Support
Size : 3,5,7,9 & 14 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



Distal Medial Humeral Locking Plate, 2.7mm (L&R)

SSL16 / TT40

Size : 3 to 9 holes & 14 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



Distal Tibia Locking Plate, 2.7mm (L&R)

SSL20 / TT44

Size : 4 to 14 holes (2 holes var.)

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



Distal Radius Plate 2.4mm

SSL39 / TT64

Type : Dorsal
Size: 5 & 6 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



L- Distal Radius Plate, 2.4mm (L&R)

SSL40 / TT65

Type : Dorsal
Size: 3+2, 4+2 holes & 3+3, 4+3 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



L- Distal Radius Plate Oblique Angled, 2.4mm (L&R)

SSL41 / TT66

Type : Dorsal
Size : 3+3 & 4+3 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



T- Distal Radius Plate, 2.4mm (L&R)

SSL42 / TT67

Type : Dorsal
Size : 3+3 & 4+3 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



Distal Radius Plate, 2.4mm (L&R)

SSL43 / TT68

Type : Volar
Size: 3+5 & 5+5 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



Buttress Locking Plate, 2.4mm (L&R)

SSL44 / TT69

Type : Volar
Size: 3+5 hole

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



Distal Radius Plate Extra-Articular, 2.4mm (L&R)

SSL45 / TT70

Type : Volar
Size: 3+5, 5+5, 5+4 holes

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



Distal Radius Plate Extra Long, 2.4mm (L&R)

SSL46 / TT71

Type : Volar
Size: 4+8 hole, 4+10 hole, 4+12 hole

Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Fixation with 2.4mm locking head screws and 2.7mm cortex screws.



**Locking Head Screw, 2.4mm
(With Self Tapping)**

SSL47 / TT78

Size: 6 to 50mm (2mm var.)

Specification

▫ Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.



**Locking Head Screw, 2.7mm
(With Self Tapping)**

SSL48 / TT79

Size: 6 to 30mm (2mm var.)

Specification

▫ Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.



Corporate Office

NRV OrthoTech Pvt. Ltd.

C-506, Signature 2 Business Park,
Near Sanand Circle, Sarkhej,
Ahmedabad - 382210, Gujarat, INDIA.
Ph : +91-79-66173506
E-Mail : info@nrvortho.com
www.nrvortho.com